

yuwell

HEAVY DUTY SUCTION APPARATUS



7A-23B type of electric suction apparatus adopts completely plastic panel design, which makes it more fancy and fashionable. The device of intercepting oil for circumfluence can reduce the oil pollution. Easy to carry, low noise, high negative pressure and large flux. It can be widely applied in the surgical operations which need to absorb phlegm, thick liquid etc.

-Oil-free pump -Convenient - no need to add oil for maintenance
-Clean - no pollution of oil -Safe - no positive pressure during using

Technical data:

Power Voltage: AC220V \pm 22V, 50Hz \pm 1Hz	Max negative pressure: \geq 90kPa
Operation noise: \leq 65 dB (A)	Input power: 180VA
Pumping frequency: \geq 40L/min	Reservoir capacity: 2500mL/pc, 2pieces