

SINGLE SCRUB SINK **EMC IV**



Dimension	
L	24 inches
W	23 inches
H	34 inches
47 inches including faucet	

Sink Dimension	
L	18 inches
W	15 3/4 inches
H	9 1/2 Inches



FOOT PUMP CONTROL

THE EMC IV SINGLE SCRUB SINK IS MADE OF HEAVY GAUGE STAINLESS STEEL SHEET. MOUNTED ON STAINLESS STEEL TUBE FRAME. SUPPLIED WITH FAUCET AND A DRAIN. THE IS EQUIPPED WITH FOOT OPERATED CONTROL FOR A ONE STEP EASY USE.